




Smart bracelet

Instruction manual

Product details



10

Plateforme d'adaptation et exigences

- 1. Android 4.4 and above
- 2. iOS 8.5 and above
- 3. Support Bluetooth 4.0 and above

Introduction of the function

- 1.Dial face interface: Three dial faces are available. Local date and time, remaining battery are shown.
- 2.Data interface: record exercise time, step count, mileage, calories burnt, sleep quality and sleep time.
- 3.Heart rate interface: the length of the heart rate page in the hand ring can measure your current heart rate.The latest 10 heart rate measurements can be displayed on this page.
- 4.Blood pressure interface: the blood pressure page of the bracelet can measure your current blood pressure.The page can display the latest 10 times blood pressure measurements.
- 5.Blood Oxygen interface: long press in the blood oxygen page can measure your current blood oxygen value.This page can display the data of the last 10 times blood oxygen tests.
- 6.Training mode: on training mode interface, press and hold to enter, there are eight training modes to choose, including walking, running, riding, jumping rope, press and hold to enter to

11


- start recording exercise time and calories burnt.Pause or stop by pressing on the fitness tracker.
- 7.Message interface: Turn on notification on APP, notification will be shown on smart bracelet wherever QQ/WeChat notification, incoming, messages are received.
- 8.Music control: after connecting with device, smart bracelet is able to control the music player. Press and hold to enter the music control panel, to play/pause, previous song, and next song.
- 9.Weather interface: The bracelet displays real-time weather conditions on the dial interface. This function works must connect to the APP on the phone must open GPS positioning and allow the APP to have access to the phone's location.If you disconnect or turn off your phone's GPS for a long time, the weather information will not be updated.
- 10.Setting: on this page, long press to select other function settings, including stopwatch/brightness adjustment/shutdown.
- 11.Others: sedentary reminder, charging remind, update notification, incoming caller ID, alarm, language matching and more.

12

APP connection method



WearFit APP QR code
(WearFit allows the use of mobile app features)




1. Scan the QR code above to download or search "Wearfit" in Android store or APP store to download and install.
 2. Press and hold the touch button for 3 seconds to turn it on, open and enter the "WearFit" app, search for the device according to the APP connection wizard, select the device model, and complete the connection pairing.
- Tips: For iOS devices, please press pair the Bluetooth for the first time connecting. Afterward, the bracelet will be able to receive the income calls, texts and more notifications pushed from the iOS.

13

Wearfit 2.0 introduction

1. Home page, additional functions, personal information



- 2.Step module: A chart will be drawn to show the steps, distance and calories burnt accordingly. Provide an insight for user of daily, weekly, and monthly data.
- 3.Sleep module: Record the daily, weekly, and monthly sleep information in the form of chart. Base on the bracelet data to calculate the daily sleep quality and sleep time.
- 4.Heart rate module: Provide the detail information of your heart rate daily, weekly, and monthly. Your heart rate will be read in every hour.

14



5. Blood pressure module: Provide the detail information of your blood pressure daily, weekly, and monthly. Your blood pressure will be read in every hour.




6. Blood oxygen module: Provide the detail information of your blood oxygen daily, weekly, and monthly. Your blood oxygen will be read in every hour.



7. Fatigue module: Measure your real time fatigue and display the information for you hourly.

15

Charging method



Please follow the instruction to charge the bracelet with charging clip when battery low.

16

Basic parameters

Screen size	1.3 inch
Bluetooth	Bluetooth4.0
Waterproof level	IP67
Type of battery	Lithium polymer battery
Battery capacity	170mAh
Charging time	2 Hours
bracelet size	260mm*40.8mm*12mm
Charging	Charging clip, Voltage 5V
Touch screen	Full screen touch
Package	Smart bracelet+Strap+Charging clip+Instruction manual

17

Warning

Please follow the doctor's instructions according to the results of self-diagnosis and treatment, it is very dangerous to have blood circulation disorders, blood disease users, please under the guidance of a doctor for treatment, the measurement results of this product are reference only not for any medical purposes and basis.

18