

PPG+ECG Blood pressure & Heart rate monitor Smart bracelet



User Manual

EN



This bracelet is USB charge directly, no need USB cable, take out the belt of touch key side, and directly insert to USB devices to charge, for example computer, power bank, USB connector.

01

Wear Guidance

The best wear position is 1-2CM behind your wrist ulnar styloid, please make sure to wear the sensor close to wrist skin, in order to get accurate results.



02

Correct testing method

Put the middle finger totally touch with the side ECG metal point, and use index finger tightly clip another side, make sure the sensor close to wrist skin, and not leak green light, then keep moveless.



03

APP download, installation and Bluetooth connection

Use mobile phone to scan the following QR code, download and install, or search "Mecare" directly from Google Play Store and Apple store, then download and install.



Mecare



Android



IOS




GooglePlay

04

After installation, open APP, register and enter personal information (make sure to enter the correct personal information, in order to get accurate blood pressure data), open mobile phone Bluetooth, "click to connect - refresh - choose the right Bluetooth ID" (Each bracelet has a sole Bluetooth ID, please reference to ON/OFF page of the bracelet)

The mobile phone will receive a pair request, please agree to connection.

After finished synchronization, the bracelet will automatic sync the time, date, week from mobile phone, and APP will sync sports data from bracelet.



05

Bracelet Functions

1.ON/OFF (MODEL NO, Version, Bluetooth ID)

Bluetooth ID is use for Bluetooth connection




2. Time, Date, Week, Battery status

For the first time use, please pair with mobile phone first in order to synchronize the correct time, date, week from mobile phone.

3. Steps



4.Distance



5.Calorie



06


6.Sport mode

In this menu, hold and long press the key to enter sports mode, it will automatic change into dynamic heart rate, and you can also see the sports distance, calorie and sport time. The heart rate records will not save to APP.

During this mode, short touch the key can switch to the time menu and then automatic switch to sport mode in 1 minute. When finished sports, hold and long press the key to quit from the mode.


Attention:

When entering sport mode, make sure to wear the bracelet 1CM behind the wrist ulnar styloid.



7.Heart rate


Turn to this menu, the bracelet will automatically start to test heart rate, when finished testing, the bracelet will vibrates, and the results will not save to APP.



07

8. Blood pressure (PPG+ECG)

Turn to this menu, the bracelet will automatically start to test blood pressure, when the icon start flicker, please put the middle finger totally touch with the side ECG metal point, and use index finger tightly clip another side, make sure the sensor close to wrist skin, and not leak green light, then keep moveless. Until the bracelet vibrates and come out the results. The results will automatic save to APP.



Attention

In order to get accurate results, please pay attention to the following points before testing

1. Before testing, make sure to enter and save correct user information in APP, such as age, sex, height, weight.
2. Before testing, please sit still and keep moveless for 5-10 minutes, make sure your heart rate keeps in a smooth condition, alarm and also smooth breathing.
3. During testing, no talking and not move, it's not suitable to test after sports or walking.
4. Keep right sitting posture, keep the hand straight with the heart.

9. Refuse calls by bracelet

After the bracelet connected with mobile phone, it will synchronize calls from mobile phone, you can refuse calls by hold and long press touch key.




08

APP specifications and functions

1. Steps

Click to check history steps, distance, calories and share to social platform. Picture(1)



(1)



(2)

2. Sleep monitor

Click to check history sleep data, such as quality, duration, wake times, fall asleep time, awake time, deep sleep time, light sleep time and share to social platform. Picture(2)

09

3. Heart rate

Click to check real-time heart rate data, you can set real-time intervals from the APP settings. Picture(3)
Also you can share data to social platform.



(3)



(4)

4. Blood pressure

Click to check the manually test records from smart bracelet and share to social platform. Picture(4)

10

5. ECG detection

Turn to this menu, the bracelet will automatically start to test blood pressure, when the icon start flicker, please put the middle finger totally touch with the side ECG metal point, and use index finger tightly clip another side, make sure the sensor close to wrist skin, and not leak green light, then keep moveless until the bracelet vibrates and come out the results. The results will automatic save to APP.



Attention


In order to get accurate results, please pay attention to the following points before testing

1. Before testing, make sure to enter and save correct user information in APP, such as age, sex, height, weight.
2. Before testing, please sit still and keep moveless for 5-10 minutes, make sure your heart rate keeps in a smooth condition, alarm and also smooth breathing.
3. During testing, no talking and not move, it not suitable to test after sports or walking.
4. Keep right sitting posture, keep the hand straight with the heart.

11

6. Running

Click to start running, the bracelet will automatic start sports mode, during running, the APP will show locus, speed, pace, calorie, distance, real-time heart rate, sports time. After finished sports, click to stop, then the records will automatic save to APP.



12

Settings

User Info: Age, Sex, Height, Weight

Connected device: (click to make settings)

1. Notifications (For iPhone, click "ON/OFF" to open and then click "Notifications" to enter settings, set which applications that you want to remind on bracelet)
2. Sedentary (click to set sedentary time, reminder method, remind time range)
3. Real-time heart rate (click to set record time intervals)
4. Find device (click to find devices and vibrate remind)
5. Dangerous HR (click to open and set dangerous heart rate level, when reach the level, bracelet will vibrate remind)
6. Wrist brightness (click to open wrist brightness)
7. Wear position (click to choose right/left hand)
8. Factory settings (click to reset the bracelet)
9. Upgrade (click to upgrade firmware version, after finished, please make sure to restart APP and connect Bluetooth again)
10. Disconnection (click to disconnect bracelet)

Sport target: set target steps

Sleep target: set everyday sleep target time

Unit setting: set metric / imperial unit

Firmware reset: Choose device to reset device firmware

About us: check APP version

Login out: Login out APP

13

Specifications:

Chip: Nordic N52832

Display: 0.96" OLED

Material: stainless steel + Aluminum + TPU

Weight: 25g

Wrist belt size: 26CM

Belt type: Removable

Charge type: USB direct charge

Battery: 100mAh Polymer lithium battery

Charge time: 1 hour

Stand by: 3-5 days

Waterproof: IP67

ECG IC: TI1291

Bluetooth version: 4.0

Support OS: Android 4.4 or above

IOS 8.0 or above

14

Attention

This device adopted medical grade IC, so will be influence by multiple factors. In order to get accuracy results, please do make sure to follow below requirements.

1. Environment: Please testing in a spacious, ventilated, silence and brightness space.
2. Wear position: Please wear the device on 1-2CM behind your wrist ulnar styloid.
3. Wear requirements: Make sure ECG panel and sensor are close to wrist skin. It's not correct to wear loose or tight, otherwise will cause leak of light or oppress blood vessel.
4. Sitting posture: Please sit straight, back to chair, put wrist on desk, keep wrist straight with your heart.
5. Testing requirements: Please keep alarm and sit still for 10 minutes, deep breath 3-5 times, no talking or moving.
6. User information: This device use multiple dimensions algorithm, please enter and save correct personal information in APP before testing, such as Age, Height, Weight, Sex.

