

LED Red Light Fractional RF Beauty Equipment

TO SCIENTIFIC AND SAFE METHOD OF
RESISTANCE YEARS INVASION



WORKING PRINCIPLE

Fractional RF is a proven radiofrequency treatment that transforms skin by working below the surface and addresses a wide range of surface imperfections that result from age and environmental factors. Benefits include the smoothing of wrinkly or uneven skin, better definition of facial features around the eyes, jaw and neckline, and even the smoothing and toning of unsightly bulges, dimples and wrinkles on the face and body. Fractional RF radiofrequency energy works with your body, remodeling your collagen to help your skin get that smooth feel and youthful look you want. Fractional RF treats targeted problem areas – so the results are focused and effective.

TAKE CONTROL OF YOUR AGING SKIN

Aging affects us all – the way we act, the way we feel, and definitely the way we look. But with advancements in skin care, you can have a say in the way aging affects your skin – and maybe even get back to the way you looked not so long ago. T Fractional RF is a proven, unique radiofrequency treatment that can help improve the appearance of sagging or loose skin, giving you a smoother, sleeker and younger look and feel. Fractional RF works in just

one treatment with little down time and delivers a natural looking result. It's you, just younger looking and more confident.

LED & FRACTIONAL RF COMBINES:



Fractional RF is a non-invasive radiofrequency therapy that can address the look and feel of aging skin on many areas of your body. Fractional RF treatments are customized to your needs, and typically require minimal recovery time. It's a great solution if you want your skin to look sleeker, smoother and younger, but can't afford a lot of time away from work or your social life.

FRACTIONAL RF IS EFFECTIVE ON:

1.The Eyes – forehead, fine lines and even our brow line can make us look worn-down. Thermage treats the upper and lower eyelids to help you look more youthful and rested

The Face – as we age, phrases like “turkey neck”, loose jowls, sagging skin and “what happened to my jaw line?” creep in.

2.Thermage smoothes lines and wrinkles and remodels collagen for the overall health of the skin

3.The Body – aging doesn’t just affect the skin on our faces, our bodies see the effects too, with crepey, sagging skin and unwanted bulges and dimples. Thermage is a non-invasive way to smooth the skin on our bodies and even temporarily improve the appearance of cellulite

Treatment Facts and Guidelines

Duration	Fractional RF treatments are typically 45 minutes for eyes or face and up to 90 minutes for larger body areas.
Treatment Regimen	One treatment is all you will need to see and feel real results.
Results	Most patients can see and feel Fractional RF results


	<p>very soon after treatment – in as little as one day. Continued improvements and contouring can happen for up to six months as your body regenerates collagen. Your body actually produces new collagen which is how the collagen matrix thickens and creates a smoothing effect. Even better, results can last years.</p>
Comfort	<p>Be sure to ask if your doctor has the latest Thermage CPT technology – it adds a gentle vibration sensation and interspersed cooling bursts that give you even more comfort during treatment. You will feel cooling and heating sensations during each pulse. Throughout the procedure, your physician will ask for feedback on the level of heat sensation to help balance your comfort with maximum results. After treatment, your skin may be slightly red, but most people resume regular activities immediately after treatment. There's typically no recovery time or special follow-up care.</p>

Safety and Side Effects

Fractional RF is built on proven science – and safety. Thermage has the longest history in aesthetic radiofrequency treatments, over a decade and over 1 million treatments. Thermage technology is

clinically proven to be safe and effective. Even so, because all patients are different and conditions can vary, some side effects are possible.

Product diagram:



- a) The probe.
- b) LED
- c) Energy key
- d) LED red light mode
- e) ON/OFF key
- f) Hot Maggie energy bond
- g) Hot Maggie mode
- h) Power jack

How to use (steps)

1. Clean the care and maintenance area.
2. Take an appropriate amount of gel or essence and apply it to the face and
3. Press the switch button (POWER) of the skin resurface instrument to turn on the device.
4. Press the instrument LED red light mode selection button (IPL) to start the red-light steady mode.

(ADJUST) Adjust the desired mode.

5. Press the thermomaji mode selection button (RF) of the instrument to start the low-grade mode.

(ADJUST) Adjust the required energy level.

Please test the sensation on the inside of your arm before the first operation.

(Choose the right gear according to the level of stimulation and heat you can accept)

6. For different parts of the care, use the corresponding care techniques to care

After the treatment is completed, press the switch (ON/OFF) to turn off the skin resurface.

Safety precautions

1. This instrument is a special beauty care instrument for healthy skin. Please do not use it for beauty treatment.
2. Please abide by the frequency of use and the contents of the operation manual. For the first time use or sensitive skin please realize a short usage time.
3. When the user's body appears abnormal, please stop using it.

Precautions before using:

If the operation range within the last six months injection of hyaluronic acid collagen, etc., is not recommended

If the operation area has been injected with artificial fillers (such as silicone), please avoid the area

If you have cardiac disease, device, or rhythm adjustment

No metal objects are allowed to be worn on the body during the operation

Please do not wear makeup on the day of the operation or clean your face before the operation.

Matters needing attention after operation:

1. Please strengthen moisturizing and sun protection.
2. It is normal for the skin to be dry after the treatment. Just pay attention to replenishing water.
3. It is recommended to use the mask at least a few times a week.

Cleaning and maintenance.

1. After each use, please wipe the body with a cotton pad or wet tissue, do not wash with water or blister cleansing.
2. Clean the beauty head with a cotton pad dipped in pure water or special alcohol (0% ethylene propylene), balance the wipe beauty head. If you use gel, please clean the beauty head thoroughly.
3. Do not use Sharp objects or metal to clean the beauty head, so as not to damage the circuit of the dot matrix.
4. After use, place the fuselage in a cool and dry environment.