

Tactile mask with LED light AN0056 – ZM888



Instructions manual

REVIEW

Mask has always been one of the essential skin care products for beauty lovers, and it enjoys great acceptance among the public. With the development of science and technology, people are no longer limited to the general medication procedure to care for the skin. Nowadays, there are various popular forms of beauty and skin care in the world such as photorejuvenation. The international medical organization has shown that the use of specific wavelengths of light can increase the vitality of skin cells and tissues, and improve skin conditions better by approaching the mode of irradiation. Photorejuvenation is the open secret of Hollywood stars. But traditional photon skin care is widely used in large photon instruments, which is expensive. It is known as the "magic light" beyond the average consumer category. LED color light skin mask is a kind of bioactive cold light imported Korea of the latest beauty technology. It produces no heat. It uses pure low energy light and unique light modulation technology to care and treat skin. It is natural and gentle, and its healing effect is remarkable. It is low risk, no side effects and other characteristics, it is the expert in skin care, suitable for all skin types.

INTRODUCTION TO THE INSTRUMENT

LED, the new replaces the old. It uses photodynamic principle to activate the deep cells, so that the skin can metabolize better. After the light is absorbed by the skin, the light energy is converted into energy in the cell, which relaxes to strengthen the microvasculature and produces photochemical reaction enzymatic reaction with the skin, increase the activity of catalase and superoxide dismutase (SOD), which are the "sanitation" and "cleaner" of the cell, and increase the decomposition of ATP, which is the energy source of the cell.

Promote cell metabolism and synthesis, increase the content of glycogen and protein, stimulate collagen production 5 times of a large area, stimulate the secretion of basic synthetic fiber growth factor (bFGF) epidermal growth factor (EGF) in the body, make the collagen fibers and elasticity, reorganize the fibers, block and eliminate the formation of melanin, accelerate cell growth, accelerate blood circulation, stimulate the production of collagen fiber cells, increase skin elasticity, repair skin aging, acne skin, spots fade, improve firmness, relieve sunburn and other care plays an important role.

METHOD OF USE

1. Clean the face skin, let the water dry naturally, or use the absorbent towel to gently remove the water on the face skin (do not use the thick towel to wipe the face tightly, which will cause the pore coarse, or the skin become rough and lose moisture).
2. Take out the LED color light skin mask instrument, hang the left and right elastic bands, remove the power control box, connect the power line of the mask instrument to the output of the power control box.
3. Controller operating instructions: press the touch mark with your finger for 3 seconds to start the machine. The color will change each time you press. The color will automatically change to start after the eighth touch.
4. The light emitting mask is placed on the face and should be adjusted to a suitable comfortable position. Both sides pass through the ear to adjust the tension of the magic belt to fix the mask.
5. We suggest that we should do the photon beauty session as much as possible during the resting time, find a comfortable sofa or bed, lie down quietly and comfortably, close your eyes and enjoy the beauty process brought by the photon. Do not read or watch TV during the time when using the photon mask.

Use it on the sofa or bed



You CANNOT read or watch TV while wearing the mask.



6. Photon mask combined with skin care products or disposable mask: With clean face, put on the skin care product or disposable mask evenly applied, put on the photon mask and use the red light irradiation, complete the set course of treatment after the time (generally recommended between 15-30 minutes), at this time, remove the mask, remove the color mask, do the cleaning work, so as to continue to use next time.

Tips: users can choose different photon colors according to their different needs. According to your skin type, you can choose the treatment that best suits your requirements.

LED LIGHT DESCRIPTION and FUNCTION

- **Red light** The wavelength of 640nm, known as the "biological activity light". It can improve cell activity, accelerate metabolism, and promote collagen formation at the same time, with whitening and reducing skin wrinkles, increase skin elasticity, improve skin condition, repair skin damage, shrink pores, firm skin, and improve antioxidant efficacy..

- **Blue light** 465nm wavelength, can inhibit inflammation, bactericidal effect, no harm to skin conditions, can make acne inflammation significantly reduced. Have effect on acne sores, anti-inflammatory bactericidal, skin reopening, no scar and pigment.

- **Green light** 532nm wavelength, with neutralization, balance and stability, to reduce skin oil secretion, balance water and fat, effectively relieve stress, dredge lymph gland and improve swelling drainage.

- **Yellow light** After the special treatment of 580nm wavelength, it can improve the oxygen of cells, supplementary energy to skin cells, make pigment decomposition, promote lymphatic drainage, improve skin, wrinkles, skin redness treatment and improve immunity.

- **Purple light** is the frequency of red and blue light. It combines two kinds of effective photo technique, especially has a particularly good and restorative effect on the treatment and elimination of acne.

- **Cyan blue light** Can gradually improve cell energy, metabolism has a very good role in this effect.

- **White light** Penetrates deep into the skin, accelerate metabolism, diminish spots, improve fine lines and loose skin.

MATTERS REQUIRING ATTENTION

1. Avoid direct vision of the mask light source or looking at the light for a long time. Keep your eyes closed during use. If you feel the light is too strong, you can cover it with cotton or mask.
2. Recommendations with quality ingredients, natural disposable mask or acne whitening skin care products, help to absorb better.
3. Because the mask material is plastic, when you open the packaging may smell, it is a normal phenomenon, the smell will disappear soon.
4. There are LED light-emitting circuits and light-emitting bodies inside the light-emitting film, which should not be extruded, bent or pressed tightly.
5. When used in bed before going to bed, it is better to put it back into the box or at least take it out of bed after treatment.
6. When cleaning the luminous mask, it is forbidden to soak it in water, wipe it with a damp cloth, the elastic band if can be washed.
7. Do not rub the mask with corrosive chemical solvents.
8. Avoid exposure to the sun.

IMPORTANT INFORMATION

European CE certified.

In the case of using medications, consult with your doctor in case of photosensitive factors of the same or other aspects that contraindicate its use.

Do not use in cases of: pregnant or breastfeeding women, children, people with thyroid disease, people allergic to light, the use of certain medications can cause allergies. In case of illness, consult your doctor.

If you have a wound on the skin should be used in a way that avoids such wound or not used.

Applying them for longer than indicated could cause a thermal effect due to heat, creating some oxidation or photo-aging. Incorrect use could cause reddening of the skin, itching and even burns, so it is essential to follow the instructions and never extend in time or frequency.

Find the appropriate combinations between dosage, frequency and intensity of treatment to achieve optimal effects on the pathologies treated.

Start with low intensity and short times of use, with rest days between sessions. Gradually increase the variables while monitoring the effects on your skin.

Never exceed the maximum use time indicated.