

# Smart BP HR bracelet

## User manual



### 1. Requirement of phones

 Android 4.2 or above  IOS 8.0 or above  Support BT4.0

### 2. Before Use

Please make sure that the electricity of the wristband is sufficient before your first use. If the electricity is too low to boot up, please charge for it first.

### 3. Download "LEFUNG HEALTH" APP



Method 1:  
Scan the QR code (For Android users)

Method 2:  
Search "LEFUNG HEALTH" from Apple Store or  
Android Market

### 4. Bind APP

Long touch the key to boot up, please make sure the BT of your phone is open

— Open "LEFUNG HEALTH" APP in your phone — Click the menu bar on the top left corner — Click "Device name / connection status" — Choose MAC IMEI which matches with your band.

### 5. Smart Band Operation

After the wrist band is connected successfully, it will sync the time and date from your phone.

APP can sync the motion state, such as Pedometer, Sleep monitor, Heart rate monitor, Blood pressure, etc.

### 6. Functions

- Home page: Display time, date and power, etc
- Smart alarm clock: Set clock in APP, it will shake to remind in time.
- Pedometer: Count your steps in a day, clear in 24 hours automatically, you can check the history steps in APP.
- Calories: Count the consumption of calories, clear in 24 hours automatically, you can check the history in APP.
- Mileage: Count the mileage, clear in 24 hours automatically, you can check the history in APP.
- Heart rate monitor: Switch to the interface of Heart Rate Monitor, it will test automatically, and save the dates in APP.

● Blood Pressure Monitor: Switch to the interface of Blood Pressure monitor, it will test automatically, and save the dates in APP.

● Sleep Monitor: Open the function of sleep detection in Band, detect your sleeping status in 24 hours. Save the dates in APP.

● Sedentary reminder: In "Band function setting", Turn on the function of Sedentary reminder and setting the time, it will shake in time.

● Drinking remind: In "Band function setting", Turn on the function of drinking reminder and setting the time, it will shake in time.

● Smart anti-lost: In "Band function setting", Turn on the function of smart anti-lost, band will shake when the phone is over the distance of bluetooth

● Message notification: After match successfully, when the phone receive the message of QQ, wechat, SMS, the band will bright and shake.

● Incoming call remind: After match successfully, when there is incoming call, band will show the contact message to remind.

● Shake for selfies: In "Band function setting", open this function, then click "Shake for selfie" in menu bar, shake the band to take photo.

● Search device: Click "Search device" in the menu bar, band will shake to remind.

The company reserves the right to make changes to the specification without notice.